Neutron flux measurement for beam polarization study at FnPB of SNS

Hitesh Rahangdale University of Tennessee Nab Collaboration





The work presented is in part supported by U.S. D.O.E award number DE-FG02-03ER41258

Neutron beta Decay and N_{ab}

$$\begin{bmatrix} \vec{v}_{e} & \vec{v}_{e} \\ \vec{v}_{d} & \vec{v}_{e} \\ \vec{v}_{w} & \vec{v}_{w} \end{bmatrix} = \begin{bmatrix} \frac{dw}{dE_{e}d\Omega_{e}d\Omega_{\nu}} \propto p_{e}E_{e}(E_{0} - E_{e})^{2} \begin{bmatrix} 1 + a \frac{\vec{p}_{e} \cdot \vec{p}_{\nu}}{E_{e}E_{\nu}} + b \frac{\vec{p}_{e}}{E_{e}} + \langle \vec{\sigma}_{n} \rangle \cdot \begin{pmatrix} \vec{p}_{e} + B \frac{\vec{p}_{\nu}}{E_{\nu}} + \cdots \end{pmatrix} \end{bmatrix}$$
Neutrino-Electron-Correlation
$$a = \frac{1 - |\lambda|^{2}}{1 + 3|\lambda|^{2}} \qquad \lambda = \frac{g_{A}}{g_{V}} \qquad A = -2\frac{|\lambda|^{2} + \text{Re}\lambda}{1 + 3|\lambda|^{2}}$$

Neutron lifetime $\Gamma = \tau_n^{-1} = \frac{2\pi}{\hbar} G_F^2 V_{ud}^2 (1+3|\lambda|^2) \int \rho(E_e) dE_e$, Fierz interference term b = 0

Current PDG valuesNab Goal: $a = -0.1049 \pm 0.0007$ $\Delta a / a = 0.1\%$ $b = 0.017 \pm 0.020$ $\Delta b = 3 \times 10^{-3}$

Uncertainty budget

	Experimental Parameter	$(\Delta a/a)_{syst}$
	Magnetic Field	6.0 · 10 ⁻⁴
Discussed in this session This Talk	Electric Potential Inhomogeneity	5.5 · 10 ⁻⁴
	Adiabaticity of Proton Motion	1 · 10 ⁻⁴
	Detector Effects	7.1 · 10 ⁻⁴
	Electron TOF	small
	Residual Gas	3.8 ·10 ⁻⁴
	TOF in Acceleration Region	3 · 10 ⁻⁴
	Background/Accidental	< 1 · 10 ⁻⁴
	Neutron Ream	3.2 · 10 ⁻⁴
	position	$1.7 \cdot 10^{-4}$
	profile	2.5 · 10 ⁻⁴
	Doppler Effect	small
	Unwanted beam polarization	1 · 10 ⁻⁴
	Total	1.2 · 10 ⁻³

10/30/22

Effect of Neutron Polarization



10/30/22

Neutron beam polarization measurement

The capture Cross section for reaction $3He + n \rightarrow p + 3H$ highly spin dependent.It is minimum when the spins are aligned and maximum with spins anti-aligned.

Transmission of Neutrons: For unpolarized He cell. For polarized cell. aligned spin For polarized cell. opposite spin

Assuming perfect spin flip of He,

$$T_{0} = Ne^{-\kappa} \quad \kappa = nl\sigma_{0}\frac{\lambda}{\lambda_{0}}$$

$$T_{up} = (1 - P_{n} \tanh(\kappa P_{He}))e^{-\kappa} \cosh(\kappa P_{He})$$

$$T_{dn} = (1 + P_{n} \tanh(\kappa P_{He}))e^{-\kappa} \cosh(\kappa P_{He})$$

$$P_{n} = \frac{R_{dn} - R_{up}}{\sqrt{\left(R_{up} + R_{dn}\right)^{2} - 4}} \qquad R_{up} = \frac{T_{up}}{T_{0}}, R_{dn} = \frac{T_{dn}}{T_{0}}$$

Although very efficient, Polarization of He is lost when spin is flipped

$$P_n = \frac{R_{down} - R_{up}}{\sqrt{\left(R_{up} + R_{down}\right)^2 - 4}} + \frac{\kappa P_{He}}{2}\bar{\epsilon}_{AFP} + \text{higher orders}$$

10/30/22

Neutron beam polarization measurement

- Beam flux measured in the upstream monitor for normalization.
- The neutron beam passes through the polarized 3He cell.
- Spin flipper is not used in this setup.
- Beam intensity is measured in the downstream monitor.
- During experiment the He spin is flipped
- Data from all detectors is collected in the
- Multi-channel voltage ADC.



Beam Flux measurement

Upstream Monitor

- Counting type detector with ${}^{14}N + CF_4$ gas mixture
- Able to detect individual neutrons.
- Used with preamp to get single neutron event.
- Produces narrow signals. ~500ns width

Downstream Monitor:

- Thick Monitor filled with ³He.
- Does not detect individual neutrons and needs high flux.
- Produces current proportional to neutron flux.

DAQ:

- Our DAQ contains voltage ADC, with 48 channels.
- It works in triggered mode
- Triggered by signal form SNS at 60Hz.
- For each trigger, data is sampled at 400us, for ~16ms, i.e, 40 time bins



10/30/22

Beam Flux Measurement

Can't directly incorporate, current upstream monitor in DAQ



- dynamic range
- Large integration time

Beam Flux Measurement



10/30/22

DNP2022, New Orleans

Beam Flux Measurement



Summary & Outlook

- Different types of Upstream and Downstream Monitors.
- Used TA to convert pulses into the instantaneous current.
- Fluctuations in upstream monitor is higher than downstream.
- Used the setup to do the beam polarization measurement.
- Increase the efficiency of upstream monitor.
- Change the DAQ system/downstream monitor.

Thanks!!

Prof. Nadia FominProf. Stefan BaeßlerDr. Seppo PenttiläSepehr Samiei



The Nab Collaboration



NC STATE UNIVERSITY









Universität Karlsruhe (TH)

Forschungsuniversität • gegründet 1825







NIVERSITY

IRGINIA







Main project funding:





Office of Science DNP2022, New Orleans

WKU





UNIVERSITY OF South Carolina